

Product Code	Carton size	Carton Quantity	Description	Latex?	Typical Application	Mandatory Field Brand Name
TT202E	38x29x21cm	50	2-way TIVA infusion set (Blue Caps) Total Length 200cm Lumen 1 - 'Side arm' with High-Flow Non-Return valve for gravity fed infusion. Lumen 2 - Drug Line with clamp and Pressure Activated Non-return / Anti-Syphon valve for Syringe-pump drug infusion.	NO	For single drug sedation or delivery of opioid analgesic agents.	Europlaz - 2-way TIVA 2.0 m Infusion Set
TT302E	38x29x21cm	50	2-way TIVA infusion set (Blue Caps) Total Length 300cm Lumen 1 - 'Side arm' with High-Flow Non-Return valve for gravity fed infusion. Lumen 2 - Drug Line with clamp and Pressure Activated Non-return / Anti-Syphon valve for Syringe-pump drug infusion.	NO		Europlaz - 2-way TIVA 3.0 m Infusion Set
TT252FE	38x29x21cm	50	2-way TIVA infusion 2-piece Combination set (Blue Caps) Total Length 250cm Lumen 1 - 'Side arm' with Swabble Needle free Valve and Non-Return valve for gravity fed infusion. Lumen 2 - Drug Line with clamp and Pressure Activated Non-return / Anti-Syphon valve for Syringe-pump drug infusion.	NO		Europlaz - 2 piece / 2-way TIVA 2.5m Combination Set
TT223E	38x29x21cm	50	3-way TIVA infusion set (Blue Caps) Total Length 220cm Lumen 1 - 'Side arm' with High-Flow Non-Return valve for gravity fed infusion. Lumens 2 & 3 - Drug Line with clamp and Pressure Activated Non-return / Anti-Syphon valve for Syringe-pump drug infusion.	NO	For the Administration of all routine TIVA / TCI	Europlaz - 3-way TIVA 2.2m Infusion Set
TT253E	38x29x21cm	50	3-way TIVA infusion set (Blue Caps) Total Length 250cm Lumen 1 - 'Side arm' with High-Flow Non-Return valve for gravity fed infusion. Lumens 2 & 3 - Drug Line with clamp and Pressure Activated Non-return / Anti-Syphon valve for Syringe-pump drug infusion.	NO	Ideal for Ophthalmic / head and neck procedures	Europlaz - 3-way TIVA 2.5m Infusion Set
TT253FE	38x29x21cm	50	3-way TIVA infusion 2-piece Combination set (Blue Caps) Total Length 250cm Lumen 1 - 'Side arm' with Swabble Needle free Valve and Non-Return valve for gravity fed infusion. Lumens 2 & 3 - Drug Line with clamp and Pressure Activated Non-return / Anti-Syphon valve for Syringe-pump drug infusion.	NO		Europlaz - 2 piece / 3-way 2.5m TIVA Combination Set
TT303E	38x29x21cm	50	3-way TIVA infusion set (Blue Caps) Total Length 300cm Lumen 1 - 'Side arm' with High-Flow Non-Return valve for gravity fed infusion. Lumens 2 & 3 - Drug Line with clamp and Pressure Activated Non-return / Anti-Syphon valve for Syringe-pump drug infusion.	NO	Designed for neuro / cardiac / max-fac procedures	Europlaz - 3-way TIVA 3.0m Infusion Set
TT224E	38x29x21cm	25	3-way TIVA infusion set (Blue Caps) Total Length 300cm Lumen 1 - 'Side arm' with High-Flow Non-Return valve for gravity fed infusion. Lumens 2, 3 & 4 - Drug Line with clamp and Pressure Activated Non-return / Anti-Syphon valve for Syringe-pump drug infusion.	NO	Perfect for cases where muscle relaxants / inotropes / insulin etc. are given.	Europlaz - 4 Way TIVA 2.2m Infusion Set
TT304E	38x29x21cm	25	4-way TIVA infusion set (Blue Caps) Total Length 300cm Lumen 1 - 'Side arm' with High-Flow Non-Return valve for gravity fed infusion. Lumens 2, 3 & 4 - Drug Line with clamp and Pressure Activated Non-return / Anti-Syphon valve for Syringe-pump drug infusion.	NO		Europlaz - 4-way TIVA 3.0m Infusion Set